


# Mölnlycke® Laparoscopic Suction and Irrigation Set


Truly customised laparoscopic solutions



Laparoscopic Suction and Irrigation Set is designed to enable efficient and easy aspiration and/or irrigation during laparoscopic procedures. It helps keep the abdominal cavity clear of liquids and blood, ensuring a safe and efficient surgical intervention.

- Ergonomic design. The handpiece can be operated with either hand
- Clear handpiece for unobstructed visualisation
- One-way valve prevents saline from refluxing
- Soft and low torque DEHP-free tubing with easy-peel separation
- Available in two lengths – 33cm and 45cm
- Double spike tubing for highest flexibility
- Colour and letter-coded buttons for clear identification of functions
- Proximal irrigation allows blockage flushing
- Compatible with the 10mm suction and irrigation probe
- Matte shaft eliminates cold light glare
- Distal lateral holes and an atraumatic tip to prevent damages to body tissues

Available both as sterile single-packed product and for inclusion in Mölnlycke® procedure trays.



Suction and irrigation



Suction and irrigation - handpiece

## Ordering information

Art. No.	Size cm	Size	Length	Size inch
899107	Suction and Irrigation Set Dual Spike	5mm	33cm	10
899108	Suction and Irrigation Set Dual Spike	5mm	45cm	10
899109	Suction and Irrigation Probe	5mm	33cm	10

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00  
 The Mölnlycke trademark, name and logo is registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
 © 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0246